# CORIOLIS® RECON

## PORTABLE AIR SAMPLER FOR AIRBORNE PATHOGENS DETECTION



COPIOLIS Microbial Air Sampler

- S High air flow rate for collection of all pathogens  $> 0.5 \ \mu m$
- Liquid sample, compatible with any type of analysis
- Ruggedized equipment for harsh environmental conditions
- Portable, easy-to-deploy and to set-up





### Collection of airborne biological agents, virus, bacteria, toxins and spores.

#### Validated collection efficiency by third party testing.



#### Coriolis<sup>®</sup> RECON – a ruggedized bio-aerosol sampler

Coriolis® RECON is a ruggedized bio-aerosol sampler, dedicated to CBRN Recon teams or first responders, with quick deployment in case of bio-threat suspicion. Coriolis® RECON is efficient, portable and has been ruggedized for unfamiliar environment. Coriolis® RECON offers flexibility with a choice of operating modes:

- ► Autonomous sampling: triggered by the operator wearing IPE; Coriolis® RECON can be used by Recon teams or first responders or for mobile applications to rapidly obtain a sample and identify the biological threat.
- Biological sentry mode: Coriolis<sup>®</sup> RECON can be set up in a standby mode, waiting for an order of a warning system.
- ▶ Long term surveillance: Coriolis® RECON can sample air up to 6 hours (Long Time Monitoring option).


#### Long Time Monitoring Option

Long Time Monitoring Option is available for the Coriolis® RECON. this option is ideal for the monitoring of a special area like Remote Base Surveillance.

This option adjusts automatically the volume of collection liquid for long sampling up to 6 hours. The volume of liquid inside the cone is always optimal for an efficient sampling.

#### Wet cyclone - a patented technology





Collection liquid is injected into the cone while aspirating.



Particles are pulled against the wall due to centrifugal force, and separated from air to be concentrated into the liquid.

**Technical features** 

Specifications	Coriolis <sup>®</sup> RECON
Application	Surveillance of critical areas
Principle	Wet cyclone
Collected particles size	> 0.5 µm
Air flow rate	600 L/min
Collection time	5, 10, 15 min - 6 hours with option
Watertight	Yes – IP54
Liquid sample	$20 \pm 5 \text{ ml}$
Dimensions ( L x w x h )	365 x 220 x 306 mm without cane
Weight	10 kg (20 kg with case)
Autonomy	1h (collection time)
Main supply	100 – 240 V
Operating temperature	+ 5°C to + 49° C (+41° F to + 120°F)
	$0^{\circ}$ C to + 49 $^{\circ}$ C with the winter pack option

www.coriolis-airsampler.com www.bertin.fr/en coriolis@bertin.fr

**Bertin Technologies** Washington DC 155 Gibbs Street 20817 Rockville, MD, USA Phone: +1 240 428 1024

Paris area, France Parc d'Activités du Pas du Lac 10 bis, avenue Ampère 78180 Montigny-le-Bretonneux Phone: +33 1 39 30 61 60

**Singapour** 2 Kallang Avenue CT Hub #07 - 23/24

BERTIN TECHNOLOGIES - Octobre 2013 - Conception IDSG - Photos by JP. Lagiewski or BSPF

0





Air Liquide Particles