

SPECIFICATION SHEET NARCOTICS TEST KIT (NTK) DEVICE PART No: BBINTK-AM / BBINTK-COC / BBINTK-HER

Narcotics Test Kit (NTK)

The Narcotics Test Kit (NTK) is a single-use, simple to use spot test. The NTK is provided as a paper strip with dry reagent test zones applied to the strip surface. The NTK addresses the increasing numbers of suspicious Drugs of Abuse shipments encountered at entry points, clandestine laboratory manufacture, illegal use, and negative social impact with a rapid, simple, and reliable technology designed for identification of Drugs of Abuse by first responders and frontline personnel.

Three NTK tests are available, capable of detecting:

- Amphetamines
- Heroin
- Cocaine

KEY FEATURES

- + The NTK tests provide results in less than 20 seconds.
- + Portability Each NTK weighs only 0.02 grams and fits in a shirt pocket or wallet.
- + Kits contain no hazardous or corrosive liquids and no glass vials.
- + No power requirements.
- + A minimum shelf life of 18 months.



Part of BBI Group



Contact: info@bbidetection.com Go to: www.bbidetection.com



SPECIFICATION SHEET CANNABIS TEST KIT (CTK) DEVICE PART No: BBICTK

Cannabis Test Kit (CTK)

The Cannabis Test Kit (CTK) is similar to the NTK in that it is a single-use, simple to use spot test. The CTK is provided as a paper strip with dry reagent test zones applied to the strip surface. A cotton swab is also provided for suspect residue transfer to the paper strip test zone for analysis. Positive indication is rapid, if not instant, taking less than 2 seconds. If no positive colour indication has been observed after 10 seconds, the suspect residue is unlikely to contain detectable cannabis material.

The CTK is capable of detecting:

- All active cannabis plant residue (e.g., leaves, stems, fruiting bodies)
- Cannabis extracts (oils, resins, and plant extract)
- Cannabis extracts within consumables, including confectionary and food products (edibles)
- Compounds including, but not limited to, cannabidiol, cannabichromene, THC, and cannabinoid

KEY FEATURES

- + The CTK test provide results in less than 10 seconds.
- Portability Each CTK weighs only 0.02 grams and fits in a shirt pocket or wallet.
- + Kits contain no hazardous or corrosive liquids and no glass vials.
- + No power requirements.
- + A minimum shelf life of 18 months.



Distributed by BBI Detection, 8 Tom McDonald Avenue, Medipark, Dundee, DD2 1NH www.bbidetection.com Email info@bbidetection.com UK facility: T: +44 (0)1382 569 900 US facility: T: +1 (800) 423 8199 Manufactured by S2 Detection Technologies LLC, 1001 16th St., Suite B-180-59 Denver, Colorado 80264 USA. © 2016 BBI Detection Ltd. All rights reserved. Patent Pending ART NTK/CTK 01

BBI Detection